



medical  
systems

# °B Connected

24/7 real-time monitoring solution  
for the safety of samples



- 24/7 real-time monitoring of internal temperature, door status, alarms, and warnings
- Monitor temperatures ranging from  $-99^{\circ}\text{C}$  to  $+99^{\circ}\text{C}$ <sup>1</sup>
- Web-based portal.<sup>2</sup> Compatible with web-based and mobile browsers
- Instant alerts via SMS and email
- Easy data access and export in the form of .pdf, .csv, .xml and .xlsx reports
- Raw data access via REST API
- Unlimited data storage and history<sup>3</sup>
- Audit trail functionality
- Also compatible with non-B Medical Systems devices<sup>4</sup>
- Available in English, French, Spanish, German, Russian and Czech

Developed according to  
GAMP® 5

1. Based on the sensors  
2. Installed on a server provided by the customer  
3. Depends on the data storage capacity of the customer's server  
4. Via temperature sensors connected to a Data Communication Unit. Conditions apply

**SAVING LIVES** THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems is a global manufacturer of medical grade refrigeration and transport solutions with over **40 years of experience**.



**B Medical Systems S.à r.l**

17, op der Hei, L-9809, Hosingen, Luxembourg



[www.bmedicalsistemas.com](http://www.bmedicalsistemas.com)

Tel: +352 92 07 31 1 | Email: [info@bmedicalsistemas.com](mailto:info@bmedicalsistemas.com) | [cit@bmedicalsistemas.com](mailto:cit@bmedicalsistemas.com)

FL\_LUX\_MD\_BC\_EN\_ED062025